



University of St Andrews



# Kinetic and mechanistic studies of the hydroformylation of 1-octene with $[\text{RhH}(\text{PPh}_3)_3\text{CO}]$ as a catalyst.

Lynzi M. Robb and David J. Cole-Hamilton

School of Chemistry, The University of St. Andrews, Purdie Building, North Haugh, St. Andrews Fife, KY16 9ST



University of St Andrews

The aim of this poster is to show that reaction progress kinetic analysis (RPKA) <sup>[1]</sup> can be modified to elucidate kinetic information from, catalysed reactions that incorporate gases.

- Gas uptake signals have been utilised as effective *in-situ* methods in RPKA. <sup>[2]</sup>
- Reactions incorporating more than one gas have not been studied in this way.
- RPKA can give a complete kinetic and hence mechanistic picture of a reaction from a minimal number of experiments.
- The reaction focused on is the hydroformylation of 1-octene using  $[\text{RhH}(\text{PPh}_3)_3\text{CO}]$ .

## Kinetic study of the hydroformylation of 1-octene using $[\text{Rh}(\text{CO})_2(\text{acac})]$ and $\text{PPh}_3$

A known reaction was studied to determine whether gas uptake can be used as a principle technique for investigation. The reaction chosen was the hydroformylation of 1-octene using  $[\text{RhH}(\text{PPh}_3)_3\text{CO}]$  because:

- Well established and accepted dissociative reaction mechanism, (fig 1).<sup>[3]</sup>
- Reaction contains three simultaneously changing reacting components.
- Results could be compared to previous kinetic and mechanistic studies in the literature.

A homogeneous batch reactor was used throughout this study (fig 2):

- Ability for efficient gas/liquid mixing.
- Can run both constant pressure (open) and constant volume (closed) experiments.
- Fitted with PicoLog software to allow for continual and accurate monitoring of pressure.





Figure 1: Dissociative mechanism for the hydroformylation of alkenes



The data collected is exported into a curve fitting software package, Igor (fig 3):

- Easy way to manipulate and present data
- Allows data to be fitted to several mathematical models

Figure 2: (left) Homogeneous batch reactor setup capable of efficient gas/liquid mixing and the ability for kinetic measurements at constant pressure and volume.  
Figure 3: (right) Screen shot of Igor showing a: the imported data (in table one) b: data being manipulated (text box, bottom left) c: data presented in graphs (graphs 3 and 0) and d: two examples of data being fitted to mathematical models; graph 3 shows and exponential fit and graph 0 shows a linear fit.



### Order in 1-octene and total gas pressure

- Open and closed experiments ran to determine the order in substrate and overall gas pressure




Figure 4: Graph showing the change in substrate concentration as a function of time for the hydroformylation of 1-octene at 40 bar and 80 °C in a closed system.

- Open system shows first order in substrate concentration




Figure 5: Graph showing the natural log of the change in substrate concentration as a function of time for the hydroformylation of 1-octene at 40 bar and 80 °C in a closed system.

- Closed system shows zero order in overall gas pressure
- Previous studies show partial orders in both hydrogen and carbon monoxide pressure <sup>[4]</sup>

### Order in hydrogen and carbon monoxide

- Series of experiments run to determine order in gas pressure
- 5 reactions run with  $\text{H}_2$  pressure between 5 and 25 bar
- Closed system used to calculate order with respect to  $\text{PCO}$




Figure 6: Graph showing the natural log plot of the rate of reaction as a function of time for the natural log of the carbon monoxide pressure  
Figure 7: Graph showing the natural log plot of the rate of reaction divided by the octene concentration and carbon monoxide pressure to their respective powers as a function of time for the natural log of the hydrogen pressure

- To check to see if opposite experiments are comparable
- 5 reactions with  $\text{CO}$  pressure between 5 and 25 bar




Figure 8: Graph showing the natural log plot of the rate of reaction as a function of time for the natural log of the hydrogen pressure  
Figure 9: Graph showing the natural log plot of the rate of reaction divided by the octene concentration and hydrogen pressure to their respective powers as a function of time for the natural log of the carbon monoxide pressure

- Results show  $\text{Rate} \propto [\text{1-octene}]^1 [\text{PH}_2]^{0.31} [\text{PCO}]^{-0.34}$  shown by 7 experiments: open, closed and 5 variable pressure experiments.
- Combining variable pressure data with closed system data allows for calculation of order in  $\text{CO}/\text{H}_2$  pressure (shown experimentally)

We gratefully acknowledge E.P.S.R.C and EastChem for funding. Peter Pogorzelec is thanked for help using the experimental apparatus throughout this work, Paul Webb is thanked for answering endless questions concerning Igor. The DCH group are thanked for all their help and support.

- 1) Blackmond, D.G., *Angewandte Chemie: International Edition*, 2005. **44**: p. 4302-4320.
- 2) Visentin, F., *et al.*, *Ind. Eng. Chem. Res.*, 2006. **45**: p. 4544-4553.
- 3) Evans, D., *et al.*, *J. Chem. Soc. (A)*, 1968, p.3133
- 4) Kamer, P.C.J., *et al.*, *Coord. Chem. Rev.*, 2004. **248**, p. 2409-2424